




## *nanoPSD*



# REAL-TIME DIGITAL SPECTROMETER WITH AN EMBEDDED PULSE-SHAPE ANALYZER AND PULSE- SHAPE DISCRIMINATOR

Model Number: PS1000

## I. FEATURES

- Real Time Operation
- Built-in PMT preamplifier with selectable sensitivity (8 settings)
- Time-Invariant Signal Processing with Zero Dead Time
- Memory Acquisition (the only system dead time ) 64ns/pulse (>15mln records/s)
- Incoming counting rates - detector signal-generation limited (> 5mln cps, LaBr)
- Time-Invariant Pulse-shape Signatures (TIPS) spectra
- One TIPS and 3 Pulse Height Spectra with 4096 channels each
- ROI selection of TIPS regions for PSD
- Simultaneous Amplitude and Pulse-Shape Discrimination

- Enhanced Pile-Up rejection and threshold settings
- Timers with typical accuracy of less than 10ppm (optional 0.2 ppm)
- USB Powered, Power Consumption 900mW (typ)
- Exceptional Temperature Stability: Gain  $< \pm 12$  ppm/°C, Base Line  $< 1$  ppm/°C.
- Temperature Operating Range: -10°C to +60°C.
- Weight <130g.
- Dimensions 3.6" x 1.5" x 1" (92 mm x 38 mm x 25 mm).
- Free *labZY-PSD* software for configuration, PSD settings, spectra acquisition and basic analysis.

## II. DESCRIPTION

The nanoPSD is part of the nanoMCA family of high-performance multichannel analyzers and radiation spectrometers. nanoPSD is real-time digital spectrometer with an embedded pulse-shape analyzer (PSA). The spectrometer can be used with scintillation detectors coupled to a photomultiplier tube (PMT). The nanoPSD has built-in PMT preamplifier and advanced digital pulse processor which operates in real time. The PSA is based on time-invariant signal processing offering high counting rates and excellent linearity over a wide dynamic range of signals. The PSA produces Time-Invariant Pulse-shape Signatures (TIPS) spectra which provide the basis for pulse-shape discrimination. Four spectra are acquired simultaneously - one TIPS and three Pulse Height Spectra with 4096 channels each. nanoPSD is a perfect match for PSD detectors such as stilbene, plastic scintillators, liquid scintillators, phosphich assemblies, and detectors with PMT anodes connected together (T-Phosphich). nanoPSD can also be used with scintillation detectors without PSD capabilities offering counting rates in excess of few million counts per second - e.g. LaBr.

As all other labZY devices nanoPSD is fully customizable allowing optimization of algorithms specifically tailored to customer requirements. Modifying functionality and signal processing algorithms of nanoPSD is as simple as a mouse click.

### III. BLOCK DIAGRAM




Fig. 1 Functional Block Diagram of the *nanoPSD*.






Fig. 2 Preamplifier sensitivity selector a). The selector is under a small cover b) on the bottom side of the enclosure.

## IV. CONNECTIONS



a)



b)

Fig. 3 nanoPSD connections.

## V. SPECIFICATIONS

### Input E:

*Signals from PMT anode: AC or DC coupled*

*Charge Sensitivity:*

*8 position sensitivity selector (SEL = 0 to 7).*

*SENSITIVITY = (8 - SEL)\*3.2 fC/channel ±5% @ gain of 1.00 and  $2^{12}$  channels.*

*Charge Sensitivity at Gain >1: Charge Sensitivity @ Gain=1 divided by the gain.*

*Fine Gain: 1.00 to 1.20 in 65536 steps.*

*Maximum Input Offset Current: ±10µA.*

*Absolute Maximum Signal Voltage: ±5V.*

*Preamplifier Time Constant: 520ns ±20%.*

### Input D:

*Type: Digital Input, 3.3V CMOS or Analog Input 0 to +2.5V.*

*Function: Analog Input to a slow 12-bit ADC.*

*Important: Leave this input unconnected when not used. **Never apply pulse or high frequency signals to this input!***

### Input/Output R:

*Type: Digital Input, 3.3V CMOS or Digital Output, 3.3V CMOS, Open Drain or Tristate.*

*Primary Input Function: Coincidence Logic Signal, 3.3V CMOS.*

*Primary Output Function: Acquisition Synchronization between Multiple Devices.*

*Custom Function: Per customer requirements.*

*Output Drive: DISABLED, PUSH-PULL, OPEN DRAIN; STRAIGHT or INVERTED.*

### Output S:

*Type: Digital Output, 3.3V CMOS, Open Drain or Tristate.*

*Primary Function:* SCA (ROI) Counting Signal, Discriminator Signals.

*Default Output Driver:* 3.3V CMOS.

*Custom Function:* Per customer requirements.

*Output Drive:* DISABLED, PUSH-PULL, OPEN DRAIN; STRAIGHT or INVERTED.

## **Output T:**

*Type:* Digital Output, 3.3V CMOS, Open Drain or Tristate.

*Primary Function:* Timing Signals, Pile-Up rejection Signal.

*Default Output Driver:* 3.3V CMOS.

*Custom Function:* Per customer requirements.

*Output Drive:* DISABLED, PUSH-PULL, OPEN DRAIN; STRAIGHT or INVERTED.

## **Digital Pulse Processor:**

*Signal Processing:* Time Invariant.

*Sampling Period:* 8ns (Frequency 125MHz).

*Quantization:* 16 bit, including offset and pile-up head room.

*Integral Nonlinearity:* 0.006% (typ), 0.018% (max) over full scale.

*Differential Nonlinearity:* <0.1% for typical high-resolution setup<sup>1</sup>.

*Peak Detection:* labZY's proprietary digital constant-fraction timing algorithm.

*Base Line Stabilizer:* Digital, Gated High-Pass Filter with Software adjustable response.

*Main Filter Digital Pulse Shape:* Trapezoidal – standard, other shapes optional.

*Main Filter Rise Time:* 16ns to 16 $\mu$ s, adjustable in increments of 8ns.

*Main Filter Flat Top:* 8ns to 2 $\mu$ s, adjustable in increments of 8ns.

*Fast Filter Digital Pulse Shape:* PSD dependent.

*Fast Filter Rise Time:* 8ns to 2  $\mu$ s, adjustable in increments of 8ns.

*Fast Filter Flat Top(Only in special FPGA designs):* 8ns to 2 $\mu$ s, adjustable in increments of 8ns.

*Digital Signal Thresholds (main and fast filters):* Automatic or manual. Adjustment in increments of one **hard size** channel.

### **Pulse-Shape Analyzer:**

*Technique:* labZY's proprietary ballistic deficit and time interval filtering algorithm.

*Output:* Time-Invariant Pulse-shape Signature(TIPS) Spectrum.

*TIPS Spectrum:* 4096 channels.

*TIPS Gain:* 1 to 128.

### **Pulse-Shape Discriminator:**

*Discrimination Technique:* ROI window selection of the TIPS peaks.

*Discrimination Windows:* 3.

*Amplitude Spectra:* 3, 4096 channels each.

*Memory Acquisition Time (the only system dead time):* 64ns, all spectra including TIPS.

### **Coincidence Circuit:**

*Coincidence Sources:* Internal timing signal and either the delayed direct logic signal at Input R or internally generated delayed logic signal (Coincidence Pulse) triggered by the edges of the logic signal at Input R.

*Modes of Operation:* Input R as coincidence/anti-coincidence window pulse; Input R edge triggered coincidence/anticoincidence pulse.

*Internal Coincidence Signal Trigger:* Selectable positive or negative edge of Input R.

*Input R Delay:* Adjustable 8ns to 32 $\mu$ s, in increments of 8ns.

*Coincidence Window:* Adjustable 8ns to 32 $\mu$ s, in increments of 8ns.

*Internal Timing Signal:* Constant Fraction Peak Detection (Peak Detect).

*Peak Detect Pulse Width:* 8ns.

*Peak Detect Delay:* Adjustable 8ns to 32 $\mu$ s, in increments of 8ns.

*Coincidence Circuit Operation Modes:* DISABLED, DIRECT, COINCIDENCE WINDOW, ANTI-COINCIDENCE WINDOW.

## **Data Acquisition:**

*Hardware Spectrum Size (**hard size**):* 4 spectra, 4096 channels each, using smart spectrum size technology. Hard size spectra are always recorded and stored in files.

*Soft Spectrum Size (Soft Size):* Instant, distortion free size conversion for display or data processing: 512, 780, 1024, 1489, 2048, 3276, 4096 channels for each spectrum. The soft size conversion does not cause destruction of the hard size spectra which allows an instant selection of any of the available soft sizes. A single acquisition allows display and/or data processing of the spectrum as any one of the soft spectrum sizes.

*Counts per Channel:* 4 bytes, 0 to 4.3 billion.

*Time Measurement:* Real and Live timers.

*Preset Time:* Real or Live.

*Timer Resolution:* 200 ns.

*Standard Timer Accuracy:*  $\pm 10$ ppm. (Includes variations due to initial tolerance, temperature and power supply voltage)

*Metrology Timer Accuracy<sup>2</sup>:*  $\pm 0.2$ ppm (TYP),  $\pm 3$ ppm all factors, including aging

*Preset Time Resolution:* 10ms.

*Maximum Preset Time:*  $43 \times 10^6$ s or 497 days.

*Dead Time Correction Technique:* Extended Paralyzable Dead Time.

*ICR Estimation:* Counting and correction for pile-up losses in either the fast channel (standard) or the main channel.

*Pile-Up Rejection:* Time between fast discriminator pulse and labZY's proprietary advanced fast discriminator pile-up detection.

*Measurement Start Time Stamp:* Start date and time UTC or LOCAL.



*Time Stamp Accuracy:* <50ms using internet NTP servers fully supported by labZY-PSD.

*Data Backup:* Battery-less. Hard Size Spectrum and All Settings.

### **Communication Interfaces:**

*Wired:* USB(also power source), Ethernet.

*Wireless:* WiFi, Bluetooth.

### **Environmental:**

*Gain Temperature Stability:* < 12 ppm/°C (typical), 20 ppm/°C (maximum)

*Base Line Temperature Stability:* Digitally stabilized, not subject to temperature drift. For comparison purposes with analog systems < 1 ppm/°C.

*Operating Temperature Range:* Normal Temperature Range -10°C to +60°C

### **Power:**

*Power Supply:* Required for all interfaces other than USB: 5V@1A wall plug or a 5V battery unit.

*Power Supply Voltage:* +5V ±10%.

*Operating Power (typ) :* 900mW at 25°C and USB interface. 800mW to 1.2W over the full Temperature Range.

*Additional Power Requirements:* nanoWF Interface - 500mW, nanoET Interface – 900mW.

Note 1: Differential Nonlinearity depends not only on the quantization properties of the digitizer, but also upon the noise level of the signal. Reference: V.T. Jordanov and K.V. Jordanova, "Quantization Effects in Radiation Spectroscopy Based on Digital Pulse Processing ", Nuclear Science, IEEE Transactions on, Vol 59, Issue 4, pp 1282 - 1288, Aug. 2012.

Note 2: Special Order.

**Mechanical:**

*Dimensions:* 3.6" x 1.5" x 1" (92 mm x 38 mm x 25 mm).

*Weight:* 135 g.




Fig. 4 nanoPSD dimensions.

## VI. APPLICATION INFORMATION

### Connecting nanoPSD to a scintillation detector:




Fig. 5 Connection diagram of the nanoPSD to a scintillation detector with a photo-multiplier tube (PMT): a) DC coupled (negative high voltage) -RECOMMENDED; b) AC coupled (positive high voltage). For optimal performance it is recommended to use a connection length of 40cm or less.

**Timing diagram of the coincidence circuit:**




Fig. 6 Timing diagrams of the built-in coincidence circuit: a) Input R as direct coincidence signal, active high or anti-coincidence signal, active low; b) Input R as direct coincidence signal, active low or anti-coincidence signal, active high.; positive edge c) and negative edge d) coincidence/anti-coincidence triggered signals.

## FPGA Design Files:

labZY provides standard FPGA designs that can be uploaded to the nanoPSD using the FPGA programming utility of the labZY-PSD software. Each version of the FPGA design comes in different files corresponding to different modes of operation of nanoPSD. Fig. 7 shows the naming specification of the FPGA design files.




Fig. 7 Naming specification of the FPGA design files.

## VIII. ORDERING INFORMATION

### **nanoPSD Pulse-Shape Analyzer and Discriminator Package PS1000**

- One **nanoPSD**, Part Number: **PS1000**

Including the following accessories:

- One USB Cable, Part Number: **NA0511**
- One BNC male to MCX male cables, Part Number: **NA0512**
- One BNC male to MCX male cables, Part Number: **NA0514**
- One Flash Drive with software and documentation

## VIII. ACCESSORIES

### **BNC female to MCX male Adapter**

Part Number NA0513

*Length: 8cm*



### **BNC male to MCX male Adapter**

Part Numbers: NA0512, NA0514

*Length: 100cm (NA0512), 40cm (NA0514)*



### **USB Data Cable ( 3ft )**

Part Number: NA0511-1

### **USB Data Cable ( 6ft )**

Part Number: NA0511-2

### **USB Data Cable ( 15ft )**

Part Number: NA0511-15


### **Bluetooth Interface Module**

Part Number: NA0520




### **Ethernet Interface Module *nanoET***

Part Number: NA0523



### **WiFi Interface Module *nanoWF***

Part Number NA0521



### ***nanoWF* Extension Cable ( 30cm )**

Part Number: NA0511-E12

### **Power Adapter**

( for use with *nanoET* and *nanoWF* )

Part Number: NA0510

*Voltage: 110/240V Current: 1A*

